

配置四：C-Flex 374热塑管版

袋体体积	管路	取件	货号1 (Renoldi膜)	货号2 (KA膜)
3D 50L	进: 100 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 100 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/8" C-Flex® 374 热塑管 + 母鲁尔接头 + 堵头	3D50LC2CP	K3D50LC2CP
3D 100L	进: 150 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 150 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/8" C-Flex® 374 热塑管 + 母鲁尔接头 + 堵头	3D100C2CP	K3D100C2CP
3D 200L	进: 150 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 150 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/8" C-Flex® 374 热塑管 + 母鲁尔接头 + 堵头	3D200C2CP	K3D200C2CP
3D 250L	进: 150 cm ID 1/2" * OD 3/4" C-Flex® 374 热塑管 + 软管堵头 出: 150 cm ID 1/2" * OD 3/4" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/8" C-Flex® 374 热塑管 + 母鲁尔接头 + 堵头	3D250C2CP	K3D250C2CP
3D 500L	进: 150 cm ID 1/2" * OD 3/4" C-Flex® 374 热塑管 + 软管堵头 出: 150 cm ID 1/2" * OD 3/4" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/8" C-Flex® 374 热塑管 + 母鲁尔接头 + 堵头	3D500C2CP	K3D500C2CP
3D 1000L	进: 150 cm ID 1/2" * OD 3/4" C-Flex® 374 热塑管 + 软管堵头 出: 150 cm ID 1/2" * OD 3/4" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/8" C-Flex® 374 热塑管 + 母鲁尔接头 + 堵头	3D10XC2CP	K3D10XC2CP