

2D 标准袋选型表 IV

配置四：C-Flex 374热塑管版

袋体体积	管路	取样	货号1(Renolit膜)	货号2(KA膜)
2D 250ml	进: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头	无	2D025C2NP	K2D025C2NP
2D 500ml	进: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头	无	2D050C2NP	K2D050C2NP
2D 1L	进: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头	无	2D01LC2NP	K2D01LC2NP
2D 1L	进: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头	10 cm ID 1/4" * OD 7/16" C-Flex® 374 + 无针取样器	2D01LC21P	K2D01LC21P
2D 1L	进: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/16" C-Flex® 374 + 母鲁尔接头 + 堵头	2D01LC2CP	K2D01LC2CP
2D 2L	进: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头	无	2D02LC2NP	K2D02LC2NP
2D 2L	进: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头	10 cm ID 1/4" * OD 7/16" C-Flex® 374 + 无针取样器	2D02LC21P	K2D02LC21P
2D 2L	进: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 1/4" * OD 7/16" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/16" C-Flex® 374 + 母鲁尔接头 + 堵头	2D02LC2CP	K2D02LC2CP
2D 5L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	无	2D05LC2NP	K2D05LC2NP
2D 5L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	10 cm ID 1/4" * OD 7/16" C-Flex® 374 + 无针取样器	2D05LC21P	K2D05LC21P
2D 5L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/16" C-Flex® 374 + 母鲁尔接头 + 堵头	2D05LC2CP	K2D05LC2CP
2D 10L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	无	2D10LC2NP	K2D10LC2NP
2D 10L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	10 cm ID 1/4" * OD 7/16" C-Flex® 374 + 无针取样器	2D10LC21P	K2D10LC21P
2D 10L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/16" C-Flex® 374 + 母鲁尔接头 + 堵头	2D10LC2CP	K2D10LC2CP
2D 20L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	无	2T20LC2NP	K2T20LC2NP
2D 20L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	10 cm ID 1/4" * OD 7/16" C-Flex® 374 + 无针取样器	2T20LC21P	K2T20LC21P
2D 20L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/16" C-Flex® 374 + 母鲁尔接头 + 堵头	2T20LC2CP	K2T20LC2CP
2D 20L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 3/4" 快装接口 (TC25) 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 3/4" 快装接口 (TC25)	10 cm ID 1/4" * OD 7/16" C-Flex® 374 + 无针取样器	2T20LC21T	K2T20LC21T
2D 35L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	无	2D35LC2NP	K2D35LC2NP
2D 35L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	10 cm ID 1/4" * OD 7/16" C-Flex® 374 + 无针取样器	2D35LC21P	K2D35LC21P
2D 35L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/16" C-Flex® 374 + 母鲁尔接头 + 堵头	2D35LC2CP	K2D35LC2CP
2D 50L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	无	2D50LC2NP	K2D50LC2NP
2D 50L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	10 cm ID 1/4" * OD 7/16" C-Flex® 374 + 无针取样器	2D50LC21P	K2D50LC21P
2D 50L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 软管堵头	40 cm ID 1/4" * OD 7/16" C-Flex® 374 + 母鲁尔接头 + 堵头	2D50LC2CP	K2D50LC2CP
2D 50L	进: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 3/8" HB ECS 母接头 出: 50 cm ID 3/8" * OD 5/8" C-Flex® 374 热塑管 + 3/8" HB ECS 公接头	10 cm ID 1/4" * OD 7/16" C-Flex® 374 + 无针取样器	2D50LC21K	K2D50LC21K